

SonoScape

X5

「LIGHTEN
YOUR IMAGING」



ECHOVITALIS

Information / Prix / Démonstration

+33(0)7 63 63 05 65

echovitalis@gmail.com



SonoScape's laptop styled color Doppler ultrasound system, the X5, will revolutionize the ultrasound market. Designed as the go-to ultrasound in any setting at the bedside or abroad such as in an ambulance thanks to its light

it the ideal system for use inside hospitals too. Be it the OR, emergency room or ICU, the X5 can be found and ready to use at anytime.

Efficient Experience

The new platform provides users a more efficient scanning process, due to its advanced technology. Since X5 aims to achieve an excellent imaging quality, these technologies contribute to produce images with high resolution and penetration, presenting a large amount of information.



Multi-beam processing technology



Spatial compound imaging



μ-Scan image processing technology



Pure inversion harmonic imaging




C-field beam




Advanced Technology

- 15.6 inch monitor in high definition
- Wide range folding angle and viewing angle
- Extended connector for up to 3 probes
- Noiseless system





- Quick boot up time
- Built-in battery supports 180 mins
- Integrated trolley with height adjustment
- Wi-Fi and Bluetooth wireless connection





Extensive Applications

various applications.



User Friendly

-  One button automatic optimization
-  Large image display region
-  Anti-reflective screen and auto brightness adjustment
-  User-defined operations